

We Speak Shoe

The New Razor–Thin High Heel Fashion Orthotic

The advanced Carboplast® graphite composite shell is just 2 mm thick. Includes a gorgeous and durable Ultrasuede® top cover. Suitable for heel heights from 2" to 5". The unique straight distal edge of the device helps stabilize the first and fifth rays, preventing inversion as the rearfoot is elevated above the transverse plane by the shoe. The orthotic is completely custom fabricated to your prescription.

Finally, your patients can have functionality, stability, comfort and fashion in a single device. Call for your free sample orthotic today!



Don't Forget Our Guarantee

We unconditionally guarantee, for a period of six months, that you and your patients will be completely satisfied with their orthotic devices.

If, for any reason, you or your patients are not satisfied, we will adjust or replace the devices at no cost to you or your patients. You won't even be charged for shipping!

866.678.4652

www.integrityortho.com

Integrity ORTHOTICS
honest people. honest answers.